



# Zirlux Zirconia *Cementation and Bonding Reference*

Zirlux® restorations may be cemented or bonded using conventional cements and luting materials. Resin modified glass ionomers and resin cements are recommended over zinc oxide phosphate cements to optimize aesthetic results. Inlays, onlays, and veneers must be adhesively cemented (bonded).

Ceramir® Crown & Bridge

**Doxa**

RelyX™ Luting Cement

**3M ESPE**

GC Fuji Plus®

**GC**

Panavia™ F2.0

**kuraray**

Multilink® Automix

**ivoclar  
vivadent**  
passion vision innovation



Refer to the respective manufacturer's instructions for use and proper application of any cementation/bonding materials.

