

# Zirlux 3D Print Resins

Zirlux 3D Print Resins are designed to seamlessly integrate into your existing digital workflow. We've validated these resins on leading 3D printers to save you time and ensure consistent, high-quality results.

Your workflow deserves premium materials as reliable and versatile as the printers you trust.



|  | Validation complete | Validation in process | Zirlux Model Print Resin | Zirlux Gingiva Mask Print Resin | Zirlux Surgical Guides Print Resin | Zirlux Nightguard & Splint Print Resin - Flexible | Zirlux Nightguard & Splint Print Resin - Rigid |
|--|---------------------|-----------------------|--------------------------|---------------------------------|------------------------------------|---|--|
| Asiga® MAX (385 nm)                                  | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Asiga® PRO 4k (385 nm)                               | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Ackuretta FreeShape (405nm)                          | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Ackuretta Dentiq (405nm)                             | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Ackuretta SOL Mono (405nm)                           | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Ackuretta SOL Plus (405nm)                           | ✓                   |                       | ✓                        |                                 | ✓                                  | ✓   | ✓  |
| Carbon® M1/M2 Series (385 nm)                        |                     |                       |                          | ✓                               | ✓                                  |   | ✓  |
| Carbon® M3 Series (385 nm)                           |                     |                       |                          |                                 | ✓                                  |   | ✓  |
| Desktop Health Envision One (385 nm)                 |                     |                       |                          |                                 | ✓                                  | ✓   |  |
| Desktop Health Einstein (385nm)                      |                     |                       |                          | ✓                               |                                    | ✓   | ✓  |
| Dekema trix™ Print (385 nm)                          | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   |  |
| HeyGears A2D (385 nm)                                |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| HeyGears A2D 4K (385 nm)                             |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| HeyGears ChairSide Pro (385 nm)                      |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| HeyGears Ultracraft A3D (385 nm)                     |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| HeyGears ChairSide (385 nm)                          |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| Rapid Shape D20+, D30+, D20II, D30II, D40II (385 nm) | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Rapid Shape D70, D90, D100 (385 nm)                  | ✓                   |                       | ✓                        |                                 |                                    |   |  |
| Sisma Everes UNO (405 nm)                            | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Microlay Eve (405 nm)                                | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Microlay Versus (385 nm)                             | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| W2P SolFlex 135/150/163170/350/363/650 (385 nm)      |                     |                       |                          | ✓                               | ✓                                  | ✓   | ✓  |
| W2P SolLab 335/370/450 (385 nm)                      |                     |                       |                          | ✓                               | ✓                                  | ✓   | ✓  |
| MiiCraft 125 (405 nm)                                | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   |  |
| MiiCraft Ultra (405 nm)                              | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   |  |
| SprintRay Pro (405 nm)                               |                     |                       |                          | ✓                               | ✓                                  | ✓   | ✓  |
| MoonRay® S100/D75 (385 nm)                           |                     |                       |                          |                                 | ✓                                  | ✓   |  |
| Flashforge Hunter (405 nm)                           |                     |                       |                          | ✓                               | ✓                                  | ✓   |  |
| Stratasys Origin One Dental (385 nm)                 |                     |                       |                          | ✓                               | ✓                                  | ✓   | ✓  |
| Phrozen Sonic XL 4k (405 nm)                         | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Phrozen Sonic Mini 4k (405 nm)                       | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| SHINING 3D® AccuFab-D1/D1s (405nm)                   | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| SHINING 3D® Accufab-L4D (405nm)                      | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Nexa3D NXD200 (405 nm)                               | ✓                   |                       | ✓                        |                                 | ✓                                  | ✓   | ✓  |
| Nexa3D XiP (405 nm)                                  | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Uniz Slash 2 (405 nm)                                |                     | □                     |                          |                                 |                                    | ✓   | ✓  |
| EPAX DX10 Pro (405 nm)                               |                     |                       |                          |                                 | ✓                                  | ✓   | ✓  |
| Planmeca Creo C5 (405nm)                             | ✓                   |                       | ✓                        |                                 | ✓                                  | ✓   | ✓  |
| LuxCreo Lux3+ (405nm)                                |                     |                       |                          |                                 |                                    | ✓   | ✓  |
| LuxCreo iLux Pro (405nm)                             |                     |                       |                          |                                 | ✓                                  | ✓   | ✓  |
| NewPro 3D NP1 (405 nm)                               |                     |                       |                          |                                 |                                    | ✓   |  |
| RayShape® Shape 1 (405nm)                            | ✓                   |                       | ✓                        | □                               | ✓                                  | ✓   | ✓  |
| Zortrax Inkspire 2 Clinical (405nm)                  | ✓                   |                       | ✓                        |                                 | ✓                                  | ✓   | □  |
| B9 Creation Dent XL (385nm)                          |                     |                       |                          | ✓                               | ✓                                  |   |  |
| Elegoo Mars 3 Pro 4k (385nm)                         | ✓                   |                       | ✓                        | ✓                               | ✓                                  | ✓   | ✓  |
| Anycubic Photon Mono X 6Ks (405 nm)                  | ✓                   |                       | ✓                        | □                               | □                                  | ✓   | ✓  |
| Anycubic Photon Mono M5 (405 nm)                     | ✓                   |                       | ✓                        | □                               | ✓                                  | ✓   | ✓  |
| ProgaPrint PR5 (385 nm)                              |                     |                       |                          |                                 |                                    | ✓   |  |
| Structo ST-01 (TBD)                                  | ✓                   |                       | ✓                        |                                 |                                    | ✓   |  |
| Roland® 3DXPRINT (TBD)                               | ✓                   |                       | ✓                        |                                 | ✓                                  | ✓   | ✓  |
| Original Prusa Medical One (405nm)                   |                     |                       |                          |                                 | ✓                                  | ✓   | ✓  |



Ready to optimize your workflow?

Visit [www.zirlux.com/printresins](http://www.zirlux.com/printresins) or scan QR code to learn more.