


Zirlux[®] Anterior Multi



- 35% Dentin Layer
- 30% Gradational Layer
- 35% Enamel Layer

Multi-Layered Zirconia

with more indications and strength than Lithium Disilicate.



(800) 496-9500
www.zahndental.com

Discover why Zirlux is *right for you.*

www.zirlux.com



Zirlux[®] Anterior Multi

Zirlux Anterior Multi is the latest development in dental zirconia with fully gradient shading that matches the Incisal, Dentin, and Gingiva of dentition—making it suitable for anterior restorations. The transition of layers mimics the translucency and esthetics of a natural tooth without adding the necessary steps to the workflow.

Zirlux Anterior Multi can be characterized after sintering for cases that require additional customization.

- Axially pressed zirconia followed by cold isostatic pressing using even pressure from every direction creates consistency and ensures the finest quality.
- Each disc is quality checked before final packaging to ensure the material is of the highest quality.
- Compatible with Roland, Amann Girrbach, and industry-leading 98.5mm CAD/CAM Milling Machines.
- Can be temporarily cemented using eugenol-free cement.

Indications:


- Anterior Crowns
- Posterior Crowns
- 2-3 Unit Bridges
- Telescopic Crowns
- Veneers

Physical Characteristics

Density before sintering	≥6.0 g/cm ³
CTE after sintering (25-500°C)	(10.5 ± 1.0) × 10 ⁻⁶ K ⁻¹
3-point flexural strength after sintering	600MPa
Accelerated aging surface monoclinic phase content	<5%
Chemical solubility after sintering	<100µg/cm ²
Cytotoxicity test	0 level
Radioactivity after sintering	<0.1Bq/g
Sintering Temperature	1430~1470°C (1450°C recommended)

Precise and accurate
results from CAD/CAM
production.

Available in 98.5mm shaded discs.



Natural transition of layers mimic Incisal, Dentin, and Gingiva.

Shade	98.5 x 10 mm	98.5 x 12 mm	98.5 x 14 mm	98.5 x 16 mm	98.5 x 18 mm	98.5 x 20 mm	98.5 x 22 mm	98.5 x 25 mm
A1	901-2307	901-2308	901-2309	901-2310	901-2311	901-2312	901-2313	901-2314
A2	901-2315	901-2316	901-2317	901-2318	901-2319	901-2320	901-2321	901-2322
A3	901-2323	901-2326	901-2327	901-2329	901-2330	901-2331	901-2332	901-2333
A3.5	901-2334	901-2335	901-2336	901-2337	901-2338	901-2339	901-2340	901-2341
B1	901-2342	901-2344	901-2345	901-2347	901-2348	901-2351	901-2353	901-2354
B2	901-2355	901-2356	901-2357	901-2358	901-2361	901-2362	901-2363	901-2364
C2	901-2365	901-2366	901-2367	901-2368	901-2369	901-2370	901-2371	901-2373
	\$106.99	\$115.99	\$131.99	\$154.99	\$171.99	\$177.99	\$194.99	\$201.99