

ZIRLUX® GINGIVA MASK

Mimic Nature with Ease



0.5Kg Bottle
(921-0032)

KEY FEATURES:

Simulates Gum Tissue: Accurately mimics the look and feel of natural gum tissue

Easily Trimmable: Allows for precise customization to achieve proper emergence profiles and pontic shapes

Flexible Yet Firm: Balances flexibility with structural integrity for effective use

Low Shrinkage: Maintains dimensional accuracy during and post printing

Optimized Workflow: Validated for leading 3D printers for seamless integration

Property	Procedure	Zirlux® Gingiva Mask
Tensile Strength (MPa)	ASTM D638	2
Young's Modulus (MPa)	ASTM D638	4.6
Elongation at Break (%)	ASTM D638	41%

Zirlux® 3D Print Resins



Explore the rest of
the Zirlux® 3D Print
Resins today

Discover why Zirlux® is *right for you.*
www.zirlux.com

© 2025 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. 24ZS8614

Exclusively Available from
ZAHN™
A DIVISION OF HENRY SCHEIN®